

WFI FILTER μ P REMOTE RESET

J. Alonso

2-May-2004

An ADAM-6060 module was mounted over the WFI μ P chassis.
The relay DO-0 of the module connects the WDI pin of the MAX-690 reset & power manager chip to ground when energized generating a clean reset signal to the μ P and associated peripherals.

The ADAM-6060 was configured as follow:

IP	134.171.122.203
SM	255.255.255.000
GW	134.171.122.002

To reset the WFI access the ADAM using the installed utility software and set DO-0 output for at least 4 seconds. Don't forget to reset DO-0, fail to do so will leave the μ P in a periodic reset cycling state.

A copy of this sheet was added to the WFI documentation together with a copy of the MAX-690 data sheet.